

# Tutorial 1

## X Ray Diffraction

Linda CRANDALL and Maggie HUFF

June 9, 2017

Date Performed: June 7, 2017  
Instructor: Chris Pratt

### 1 Objective

To get to know the diffractometer and X'Pert Data Collector

### 2 Results and Conclusions

Increasing the width of the slits generally increases the intensity because more photons can get through. If you decrease the size of the anti-scatter slit you reduce the background signal so the peaks get smoother but also less intense. The accepted wavelengths for  $K\alpha_1$  and  $K\alpha_2$  are .1540562 and .1544398 nm, respectively.




Figure 1: Si sample with Cu target:  $K\alpha_1$  and  $K\alpha_2$  peaks

### 3 Discussion

$CuK_{\alpha}1$  is about twice the intensity of  $CuK_{\alpha}2$  because the transition for the first one

### 4 Answers to Definitions

- a. The *atomic weight of an element* is the relative weight of one of its atoms compared to C-12 with a weight of 12.0000000... , hydrogen with a weight of 1.008, to oxygen with a weight of 16.00. Atomic weight is also the average weight of all the atoms of that element as they occur in nature.
- b. The *units of atomic weight* are two-fold, with an identical numerical value. They are g/mole of atoms (or just g/mol) or amu/atom.
- c. *Percentage discrepancy* between an accepted (literature) value and an experimental value is

$$\frac{\text{experimental result} - \text{accepted result}}{\text{accepted result}}$$