

Tutorial 1

X Ray Diffraction

Linda CRANDALL and Maggie HUFF

June 26, 2017

Date Performed: June 7, 2017
Instructor: Chris Pratt

1 Objective

To get to know the diffractometer and X'Pert Data Collector

2 Results and Conclusions

Scan 4 is the best to use for 422 planes in Silicon; it uses a time per step of .5 seconds over a 1 degree range and .02 degrees step size. The large and small peaks of $CuK_{\alpha}1$ and $CuK_{\alpha}2$ do appear to be very close to a 2:1 ratio in intensity.




Figure 1: 422 plane in Si sample with Cu target: $K\alpha_1$ and $K\alpha_2$ peaks, variable step size, scan time, time per step

3 Discussion

Increasing the scan time lets in more signal, but also lets in more background. Increasing the step size makes the peaks less noisy. Using continuous scan mode instead of step mode makes the scan time shorter.

4 Answers to Multiple Choice Questions

- (i) d
- (ii) c
- (iii) c
- (iv) b